

Technical data

Lyla sterilization programs

Sterilization cycle:		S Wrapped 134		S Prion 134		S Gentle 121		S Unwrapped 134	
Temperature:		134°C		134°C		121°C		134°C	
Holding time (minutes):		4'		18' 30"		20' 30"		4'	
Total cycle duration (minutes) including complete drying:		Automatic from	Typical load	Automatic from	Typical load	Automatic from	Automatic from	Typical load	
	Lyla 17:	38'–45'	42'	52'–59'	56'	61'–64'	20'–26'	23'	
	Lyla 22 (110 V):	49'–56'	51'	63'–71'	65'	69'–73'	28'–35'	30'	
	Lyla 22 (220 V):	39'–46'	43'	53'–60'	57'	62'–66'	21'–27'	24'	
Max. load Lyla 17/22:		max. 4 kg/typical 2 kg/max. porous 0.5 kg						max. 4 kg unwrapped	
Max. container load:		max. 4 kg							
Test cycles:		Bowie and Dick/Vacuum							
Item supplied with reversible rack, 3 aluminium trays, drain hose, tray holder, 8 GB USB stick, door opening tool, factory test									

Lyla



Series: Type:	Lyla 17 RIN-210	Lyla 22 RIN-210
Chamber size:	17 l	22 l
Power supply:	200–240 V AC; 50 Hz or 60 Hz; 10 A single phase	100–125 V AC; 50 Hz/60 Hz; 12 A single phase
		200–240 V AC; 50 Hz or 60 Hz; 10 A single phase
Power consumption:	2.0–2.4 kW	110 V: 1.2–1.5 kW
		220 V: 2.0–2.4 kW
Overall dimension (w × h × d):	465 × 452 × 646 mm	
Weight (empty):	40 kg	42 kg
Main/used water tanks:	4.8/4.8 l	
Working range:	7 to 10 cycles	
Usable space in chamber (w × h × d) ⁽¹⁾ :	195 × 195 × 312 mm	195 × 195 × 400 mm
Connection types:	1 USB port in the front (additional 1 in the rear as option), automatic water filling kit (optional)	

(1) Usable space is referred to standard rack configuration.

Lyla sterilizers were designed, certified and validated with the most stringent directives and standards:

N°2017/745	Medical Device Regulation	IEC 61326-1	Electromagnetic compatibility
2014/68/EU	Pressure equipment Directive	IEC 61010-1	Safety requirements sterilizer
EN 13060	Small steam sterilizer	IEC 61010-2-040	Specific requirements for steam sterilizer
IEC 61770	Electric appliances connected to the water mains	The sterilizer can be validated according to EN ISO 17665-1	

Manufacturer:

W&H Sterilization Srl.
Via Bolgara, 2
24060 Brusaporto (BG), Italy
t +39 035-66 63 000
f +39 035-50 96 988
office.sterilization@wh.com
wh.com

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

20550 AEN Rev. 002 / 02.03.2023
Subject to alterations