

Technical data

Figures and facts at a glance.

MELAG washer-disinfectors are available in the versions 400 Volt (three-phase current, DTA) and 230 Volt (alternating current, DTB). If the installation location permits connection to a 400 Volt supply, we recommend the DTA version due to its significantly shorter operating times.

Programs	MELAtherm® 10		MELAtherm® 10 Evolution	
	DTA / 400 V*	DTB / 230 V*	DTA / 400 V*	DTB / 230 V*
Quick-Program	36 min	53 min	37 min	57 min
Universal-Program	40 min	59 min	42 min	63 min
Intensive-Program	51 min	64 min	54 min	67 min
Ophthalmic-Program	42 min	59 min	43 min	63 min

* All programs plus active drying, which the user can interrupt at any time.

	MELAtherm® 10	MELAtherm® 10 Evolution
Device dimensions (W × H × D)	59,8 × 81,8 (83,6**) × 67,8 cm fits under a 60 cm worktop	
Device dimensions (W × H × D) with floor unit	59,8 × 124 × 67,8 cm	
Washing chamber (W × H × D)	46,5 × 40,5 × 44,8 cm	
Weight	80 kg (86 kg**)	
Weight with floor unit	106 kg	
Electrical connection DTA	3 N AC 400 V, 50 Hz, 3 × 16 A	
Power input DTA	9,3 kW	10,2 kW
Electrical connection DTB	AC 220 –240 V, 50 Hz, 1 × 16 A	
Power input DTB	3,3 kW	3,3 kW
Dimensions MELAdem® 53 / 53 C	Ø 24 × 57 cm (61,5 cm inc. connecting bend) / Ø 24 × 44,3 cm (48,8 cm)	

** with cover plate



Save paper and protect the environment:

To save paper, we provide the user documentation electronically in our download center.
For further information please visit: www.melag.com

MELAG
competence in hygiene